

Take see and treat to the next level with the *new* Omni™ hysteroscope

Designed to allow you to perform your best work

3 HYSTEROSCOPES IN 1

Convenience to see & treat
pathology with a single
hysteroscope

ROD LENS


Quality visualization
throughout
procedures

LONG WORKING LENGTH (200 mm)

Facilitates easier
access and treatment
for larger patients

SMALL DIAMETER

- Easier insertion
- Greater patient comfort



Two kit options to fit your practice needs

	Standard Kit (60-250-1)	Light Kit (60-250-2)
Base Scope	✓	✓
3.7 mm Diagnostic Sheath	✓	✓
5.5 mm Sheath and Outflow Channel	✓	✓
6 mm Sheath and Outflow Channel	✓	
Light Guide Adapters	✓	✓

THE OMNI HYSTEROSCOPE MYOSURE COMPATIBILITY	MYOSURE MANUAL DEVICE	MYOSURE LITE DEVICE	MYOSURE REACH DEVICE	MYOSURE XL DEVICE
5.5 mm Sheath	✓	✓	✓	
6 mm Sheath	✓	✓	✓	✓

ORDERING	
Description	Order Number
The Omni™ Hysteroscope Standard Kit	60-250-1
The Omni™ Hysteroscope Light Kit	60-250-2
Rod Lens Hysteroscope Outflow Channel	40-201
Rod Lens Hysteroscope XL Outflow Channel	50-201XL
Rod Lens Scope Seals – 10 per pack	40-902
Hysteroscope Inflow and Outflow Channel End Cap	40-904
Hysteroscope Polishing Paste (3 Tubes)	40-905
Light Post Adapters - Wolf/Storz Light Guide Adapter Set	ASY-04996
Omni Instrument Tray	60-903

Contact a sales representative to learn more.

For detailed benefit and risk information, please refer to the IFU.

IMPORTANT SAFETY INFORMATION

The Omni™ hysteroscope is intended to provide viewing of the cervical canal and uterine cavity for the purpose of performing diagnostic and operative procedures. The Omni hysteroscope is not appropriate for patients who are exhibiting acute pelvic inflammatory disease or exhibit the following conditions: inability to distend uterus, cervical stenosis, cervical/vaginal infection, uterine bleeding or menses, invasive carcinoma of the cervix, recent uterine perforation, medical contraindications, or intolerance to anesthesia.

