


Omni™ 4K

Video System

State of the art visualization
system for **GYN Surgery.**



ADVANCED.
BRILLIANT.
COMPLETE.

CAPTURE 4X THE INFORMATION WITH OMNI 4K¹

ADVANCED.

Advanced autoclaveable camerahead with 4K resolution that results in vibrant colors, posterior illumination, and even light distribution.

BRILLIANT.

Brilliant 4K resolution that allows for four times the amount of information as HD. 4K captures a wider spectrum of color during surgery.

COMPLETE.

Complete 3-in-1 system that automatically captures and stores data. The Omni 4K console is the cornerstone of the system which contains the Camera Control Unit, the LED Light Engine, and the Digital Documentation System.

Visit Hologic.com or contact your sales representative for more information about our products.

Omni™ 4K
Video System

REFERENCES

1 As compared to HD.

ADS-03009 Rev 001 ©2020 Hologic, Inc. All rights reserved. Specifications are subject to change without prior notice. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Hologic, Omni, The Science of Sure and associated logos are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries.

ORDERING INFORMATION

OMNI 4K

PRODUCT NUMBER	PRODUCT DESCRIPTION
83-10-4010	4K UHD Camera Control Unit, Hologic
83-10-4120	4K UHD Camera Head, Integrated, 1 CMOS, Hologic
83-12-4001	4K UHD Office Monitor, Hologic
83-20-5003	5 MM Light Guide, Blue
83-20-50031	Wolf Light Guide Post Adapter
83-20-50032	Acmi Light Guide Post Adapter
83-20-50033	Storz Light Guide Post Adapter

PACKAGES

PRODUCT NUMBER	PRODUCT DESCRIPTION
OVS-100	Complete Starter Package
OVS-101	Camera and CCU Package

