

Hologic now distributing

Autoclave

Intuitive by design.



SIMPLE. FAST. EFFICIENT.

SPEED AND SIMPLICITY IN STERILIZATION

SIMPLE.

An intuitive menu interface and high-end color touch display.

FAST.

A large chamber capacity with small countertop footprint and versatile cycle program options.

EFFICIENT.

Speed of sterilization allows for rapid case turnover.

ORDERING INFORMATION

PRODUCT NUMBER	ITEM DESCRIPTION
90000284	Lexa Autoclave Value
90000286	Lexa Autoclave Complete
MN-111	User Manual



CAPACITY

11" diameter, 28 liters volume

Comes with 4 large trays

Holds up to 4 levels of cassettes placed horizontally and vertically

Visit [Hologic.com](https://www.hologic.com) or contact your sales representative for more information about our products.

Autoclave